

## Product datasheet for **AP09837PU-N**

### Testosterone Sheep Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	<b>ELISA:</b> 10 µg/ml.
Host:	Sheep
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Testosterone (17)-BSA
Specificity:	This antibody reacts to Testosterone.
Formulation:	20mM Phosphate, 150mM Sodium Chloride, pH 7.2 containing 0.09% Sodium Azide as preservative State: Ig Fraction State: Liquid Ig fraction prepared by Caprylic Acid and Ammonium Sulphate precipitation procedures
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Upon receipt, store undiluted (in aliquots) at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Background:	Testosterone is the principal androgen, or male sex hormone. One of the group of compounds known as anabolic steroids, testosterone is secreted by the testes but is also synthesized in small quantities in the ovaries, cortices of the adrenal glands, and placenta, usually from cholesterol. Testosterone is necessary in the fetus for the development of male external genitalia. Testosterone also stimulates protein synthesis and accounts for the greater muscular development of the male. For many years, synthetic steroids similar to testosterone have been used by athletes with the goal of improving performance, although recent medical research has shown that these drugs may have a wide range of harmful side effects.



[View online »](#)