

Product datasheet for **AP09839PU-N**

Tetracycline (Amide) Sheep Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA: 5.0µg/ml.
Host:	Sheep
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Tetracycline Amide-BTG
Specificity:	This antibody reacts to Tetracycline (Amide).
Formulation:	20mM Phosphate, 150mM Sodium Chloride, pH 7.2 containing 0.09% Sodium Azide as preservative State: Ig Fraction State: Liquid Ig fraction prepared by Caprylic Acid and Ammonium Sulphate precipitation procedures
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Synonyms:	Criseociclina, Enterocycline, Tsiklomistsin, Abricycline, Achromycin, Agromicina, Ambramicina, Ambramycin, Biocycline



[View online »](#)