

## Product datasheet for **AP09841PU-L**

### Thiabendazole Sheep Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	<b>ELISA:</b> 2.5 µg/ml.
Host:	Sheep
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Thiabendazole (5)-BTG
Specificity:	This antibody recognizes Thiabendazole.
Formulation:	20mM Phosphate, 150mM Sodium Chloride, pH 7.2 State: Ig Fraction State: Liquid Ig fraction prepared by Caprylic Acid and Ammonium Sulphate precipitation procedures Preservative: 0.09% Sodium Azide
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)

**Background:**

Anthelmintics or anti-helminthics are a class of drugs that are effective against a range of intestinal parasitic worms (helminths). Parasitic helminths must maintain an appropriate feeding site. Nematodes and trematodes must actively ingest and move food through their digestive tracts to maintain an appropriate energy state; these together with reproductive processes require a well defined and developed neuromuscular coordination. Anthelmintic treatment is a multi-targeting system designed to interfere with the integrity of parasite cells. The pharmacologic basis of the treatment for helminths involves the targeting of neuromuscular coordination, or protective mechanisms against host immunity, which lead to starvation, paralysis, and expulsion of the parasite. The benzimidazole class of drugs were introduced in 1961 and interfere with the parasite's ATP pathway on a cellular level. They bind to a specific building block called  $\beta$ -tubulin and prevent its incorporation into certain cellular structures called microtubules, which are essential for energy metabolism.

**Synonyms:**

TBZ, Tiabendazole, Mintezol, Thiabendazole, Mintesol, Omnicid, Lombristop, Thibenzol, Mintezol, Tresaderm, Arbotect