

Product datasheet for **AP09843PU-L**

Tobramycin Sheep Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA: 1.0 µg/ml
Host:	Sheep
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Tobramycin-BSA
Specificity:	This antibody recognizes Tobramycin
Formulation:	20mM Phosphate, 150mM Sodium Chloride, pH 7.2 containing 0.09% Sodium Azide as preservative. State: Ig Fraction State: Liquid Ig fraction prepared by Caprylic Acid and Ammonium Sulphate precipitation procedures.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody (in aliquots) at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Background:	Tobramycin is an aminoglycoside antibiotic used to treat infections caused by many different bacteria. Tobramycin is usually administered by intravenous infusion, intramuscular injection, or inhalation.
Synonyms:	Tobrex, Tobracin, Tenebrimycin, Tobramitsetin, Tobramycetin, Gotabiotic, Tobradistin, Tobramaxin, Tobramicin



[View online »](#)