

Product datasheet for **AP09844PU-N**

Tolfenamic Acid Sheep Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA: 1.25 µg/ml
Host:	Sheep
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Tolfenamic Acid-BTG
Specificity:	This antibody recognizes Tolfenamic Acid.
Formulation:	20mM Phosphate, 150mM Sodium Chloride, pH 7.2 containing 0.09% Sodium Azide as preservative. State: Ig Fraction State: Liquid Ig fraction prepared by Caprylic Acid and Ammonium Sulphate precipitation procedures.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody (in aliquots) at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Background:	Tolfenamic acid (TA) is one of the class of non-steroidal anti-inflammatory drugs. It is used to treat the symptoms of migraine.
Synonyms:	Clotam, Rociclyn, Bifenac



[View online »](#)