

## Product datasheet for **AP09845PU-L**

### Toltrazuril Sheep Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	The antibody is suitable for the development of Immunoassays or Immunoaffinity purification columns. <b>ELISA:</b> 0.15 µg/ml. 10 ng/ml Toltrazuril produces 53% inhibition in a competitive ELISA, employing Toltrazuril polyclonal antibody.
Host:	Sheep
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Toltrazuril-BTG
Formulation:	20mM Phosphate, 150mM Sodium Chloride, pH 7.2 containing 0.09% Sodium Azide as preservative. State: Ig Fraction State: Liquid Ig fraction prepared by Caprylic Acid and Ammonium Sulphate precipitation procedures.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody (in aliquots) at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Synonyms:	Baycox



[View online »](#)