

Product datasheet for **AP09846PU-L**

Tramadol (3) Sheep Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA: 1.25 µg/ml
Host:	Sheep
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Desmethyl Tramadol (3)-BTG
Specificity:	This antibody recognizes Tramadol
Formulation:	20mM Phosphate, 150mM Sodium Chloride, pH 7.2 State: Ig Fraction State: Liquid Ig fraction prepared by Caprylic Acid and Ammonium Sulphate precipitation procedures. Preservative: 0.09% Sodium Azide
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody (in aliquots) at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Background:	Tramadol, (1R,2R)-2-(dimethylaminomethyl)-1-(3-methoxyphenyl)cyclohexan-1-ol, is a narcotic analgesic used for severe pain. It may be habituating.
Synonyms:	Ultram, Ralivia flashtab, Tramadolum, Ultracet, Amanda, Tramal, Ralivia ER, Tradonal odis



[View online »](#)