

## Product datasheet for **AP10442PU-N**

### ANP (NPPA) Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	This antibody reacts with 17 kDa ANF protein clived from Natriuretic peptides A.
Formulation:	0.1M Tris, 0.1M Glycine, 2% Sucrose State: Purified State: Lyophilized purified powder
Reconstitution Method:	Restore in distilled water.
Concentration:	lot specific
Purification:	Affinity chromatography on Protein A
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	natriuretic peptide A
Database Link:	<a href="#">Entrez Gene 4878 Human P01160</a>
Synonyms:	Atrial Natriuretic Peptide, Prepronatriodilatin, NPPA, ANP, PND



[View online »](#)