

## Product datasheet for **AP10448SU-N**

### LRP2 Rabbit Polyclonal Antibody

#### Product data:

|                        |   |
|------------------------|---|
| Product Type:          | Primary Antibodies  |
| Applications:          | ELISA, IHC, WB  |
| Recommended Dilution:  | <b>ELISA.</b><br><b>Western Blot.</b><br><b>Immunohistochemistry.</b>   |
| Reactivity:            | Human   |
| Host:                  | Rabbit  |
| Clonality:             | Polyclonal  |
| Immunogen:             | Synthetic peptides derived from the C-term part of the human protein.   |
| Specificity:           | This antibody reacts with 522 kDa LRP-2 protein.  |
| Formulation:           | State: Serum<br>State: Lyophilized serum<br>Preservative: None  |
| Reconstitution Method: | Restore in distilled water.   |
| Conjugation:           | Unconjugated  |
| Storage:               | Prior to reconstitution store at -20°C.<br>Following reconstitution store undiluted at 2-8°C for one month<br>or (in aliquots) at -20°C for longer.<br>Avoid repeated freezing and thawing. |
| Stability:             | Shelf life: one year from despatch.   |
| Gene Name:             | LDL receptor related protein 2  |
| Database Link:         | <a href="#">Entrez Gene 4036 Human P98164</a>   |
| Synonyms:              | Megalin, Glycoprotein 330, LRP-2  |



[View online »](#)