

Product datasheet for **AP13802PU-N**

PANK1 (N-term) Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	ELISA: 1/1,000. Western Blot.
Reactivity:	Human
Host:	Rabbit
Isotype:	Ig
Clonality:	Polyclonal
Immunogen:	This antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide selected from the N-terminal region of human PANK1.
Specificity:	This antibody reacts to PANK1.
Formulation:	PBS with 0.09% (W/V) sodium azide State: Purified State: Liquid purified Ig
Concentration:	lot specific
Purification:	Protein G column, eluted with high and low pH buffers and neutralized immediately, followed by dialysis against PBS
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	pantothenate kinase 1
Database Link:	Entrez Gene 53354 Human Q8TE04
Background:	PANK1 belongs to the pantothenate kinase family. Pantothenate kinase is a key regulatory enzyme in the biosynthesis of coenzyme A (CoA) in bacteria and mammalian cells. It catalyzes the first committed step in the universal biosynthetic pathway leading to CoA and is itself subject to regulation through feedback inhibition by CoA.



[View online »](#)

Synonyms: Pantothenate kinase 1, Pantothenic acid kinase 1, PANK1a, PANK1b, MGC24596

Product images:



(LEFT) Western blot analysis of PANK1 antibody (N-term) in MCF-7 cell line lysates (35ug/lane). PANK1 (arrow) was detected using the purified Pab. (RIGHT) Western blot analysis of PANK1 antibody (N-term) in mouse spleen tissue lysates (35ug/lane). PANK1 (arrow) was detected using the purified Pab.