

Product datasheet for **AP14950PU-N**

PI 3 Kinase p85 beta (PIK3R2) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	ELISA: 1/1,000. Western blotting: 1/50 - 1/100. Immunohistochemistry: 1/10 - 1/50.
Reactivity:	Human
Host:	Rabbit
Isotype:	Ig
Clonality:	Polyclonal
Immunogen:	This antibody is generated from rabbits immunized with a peptide corresponding to amino acid residues surrounding Y467 of human PIK3R2.
Specificity:	This antibody reacts to PIK3R2.
Formulation:	PBS with 0.09% (W/V) sodium azide State: Aff - Purified State: Liquid purified Ig
Concentration:	lot specific
Purification:	Protein A column, followed by peptide affinity purification
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	phosphoinositide-3-kinase regulatory subunit 2
Database Link:	Entrez Gene 5296 Human O00459




[View online »](#)

Background:


Protein kinases are enzymes that transfer a phosphate group from a phosphate donor, generally the γ phosphate of ATP, onto an acceptor amino acid in a substrate protein. By this basic mechanism, protein kinases mediate most of the signal transduction in eukaryotic cells, regulating cellular metabolism, transcription, cell cycle progression, cytoskeletal rearrangement and cell movement, apoptosis, and differentiation. The family has been classified in 8 major groups based on sequence comparison of their tyrosine (PTK) or serine/threonine (STK) kinase catalytic domains. PIK3R2 binds to activated Protein Tyrosine Kinases, which are phosphorylated, through its SH2 domain, and acts as an adaptor, mediating the association of the P110 catalytic unit to the plasma membrane.

Synonyms:

PI 3 Kinase p85 beta, PI3K regulatory subunit beta

Product images:

Western blot analysis of PIK3R2 (arrow) using rabbit polyclonal PIK3R2 Antibody. 293 cell lysates (2 μ g/lane) either nontransfected (Lane 1) or transiently transfected with the PIK3R2 gene (Lane 2)



Formalin-fixed and paraffin-embedded human brain tissue reacted with PIK3R2, which was peroxidase-conjugated to the secondary antibody, followed by DAB staining.