

Product datasheet for **AP17379PU-N**

FKBP14 (N-term) Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	Western blot: 1:50 - 1:100. ELISA: 1:1,000.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	KLH conjugated synthetic peptide selected from the N-terminal region of human FKBP14
Specificity:	This antibody detects FKBP14 at N-term.
Formulation:	PBS with 0.09% (W/V) sodium azide State: Liquid Ig fraction
Concentration:	lot specific
Purification:	Prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS
Conjugation:	Unconjugated
Storage:	Store the antibody at 2 - 8 °C up to one month or (in aliquots) at -20 °C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	FK506 binding protein 14
Database Link:	Entrez Gene 55033 Human Q9NWM8
Background:	PPIases accelerate the folding of proteins during protein synthesis.
Synonyms:	FKBP-14, FK506-binding protein 14, FKBP-22
Note:	Molecular weight: 24172 Da
Protein Families:	Druggable Genome



[View online »](#)

Product images:



Western blot analysis of FKBP14 Antibody (N-term) in NCI-H460 cell line lysates (35 ug/lane). FKBP14 (arrow) was detected using the purified Pab.