

OUR RESEARCH

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for AP17686PU-N

PSME2 (C-term) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	ELISA: 1/1,000. Western blotting: 1/50 - 1/100. Immunohistochemistry on Paraffin embedded sections 1/50 - 1/100
Reactivity:	Human, Mouse
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	KLH conjugated synthetic peptide selected from the C-terminal region of human PSME2
Specificity:	This antibody reacts to PSME2.
Formulation:	PBS with 0.09% (W/V) sodium azide State: Liquid purified Ig
Concentration:	lot specific
Purification:	Saturated Ammonium Sulfate (SAS) precipitation
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	proteasome activator subunit 2
Database Link:	<u>Entrez Gene 19188 MouseEntrez Gene 5721 Human</u> <u>Q9UL46</u>



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2024 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

SME2 (C-term) Rabbit Polyclonal Antibody – AP17686PU-N


Background:The 26S proteasome is a multicatalytic proteinase complex with a highly ordered structure
composed of 2 complexes, a 20S core and a 19S regulator. The 20S core is composed of 4
rings of 28 non-identical subunits; 2 rings are composed of 7 alpha subunits and 2 rings are
composed of 7 beta subunits. The 19S regulator is composed of a base, which contains 6
ATPase subunits and 2 non-ATPase subunits, and a lid, which contains up to 10 non-ATPase
subunits. Proteasomes are distributed throughout eukaryotic cells at a high concentration
and cleave peptides in an ATP/ubiquitin-dependent process in a non-lysosomal pathway. An
essential function of a modified proteasome, the immunoproteasome, is the processing of
class I MHC peptides. The immunoproteasome contains an alternate regulator, referred to as
the 11S regulator or PA28, that replaces the 19S regulator. Three subunits (alpha, beta and
gamma) of the 11S regulator have been identified. PSME2 is the beta subunit of the 11S
regulator, one of the two 11S subunits that is induced by gamma-interferon.

Synonyms: Proteasome activator, PA28b, PA28beta


Protein Pathways:

Antigen processing and presentation, Proteasome

Product images:



Western blot analysis of PSME2 antibody (C-term) in mouse lung tissue lysates (35ug/lane). PSME2 (arrow) was detected using the purified Pab.



Formalin-fixed and paraffin-embedded mouse hepatocarcinoma reacted with PSME2 Antibody (C-term), which was peroxidase-conjugated to the secondary antibody

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2024 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US