

Product datasheet for **AP20259PU-N**

NIFK Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	IF, IHC
Recommended Dilution:	Immunohistochemistry on paraffin sections: 1/50-1/200. Immunofluorescence: 1/50-1/200.
Reactivity:	Human, Mouse
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	This antibody detects endogenous levels of NIFK protein.
Formulation:	Phosphate buffered saline (PBS) with 0.05% sodium azide, approx. pH 7.2 State: Aff - Purified State: Liquid purified Ig fraction
Concentration:	1.0 mg/ml
Purification:	Affinity-chromatography using epitope-specific immunogen; purity is > 95% (by SDS-PAGE)
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	nucleolar protein interacting with the FHA domain of MKI67
Database Link:	Entrez Gene 67949 Mouse Entrez Gene 84365 Human Q9BYG3
Background:	The structural proteins for the complex metalloenzyme nitrogenase include NIFK, NIFD and NIFH. These proteins are all necessary for archaeal and bacterial nitrogen fixation. The NIFK gene encodes the b subunit of the nitrogenase molybdenum-iron (MoFe) tetramer. NIFK localizes to the nucleolus where it interacts with the fork-head associated domain of the proliferation marker protein Ki-67 in a mitosis-specific and phosphorylation-dependent manner. NIFK is widely expressed in adult tissues, suggesting other functions in addition to its interaction with Ki-67, which is only expressed in proliferating cells.
Synonyms:	MKI67IP, NOPP34



[View online »](#)

Protein Families: Transcription Factors

Product images:



Immunohistochemistry analyzes of NFKB1 antibody (AP20259PU-N) in paraffin-embedded human lung carcinoma tissue.