

Product datasheet for **AP20400PU-M**

BubR1 (BUB1B) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	Western blot: 1/500-1/1000.
Reactivity:	Human, Mouse
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	This antibody detects endogenous levels of ubR1 protein. (region surrounding Thr368)
Formulation:	Phosphate buffered saline (PBS), pH 7.2. State: Aff - Purified State: Liquid purified Ig fraction Preservative: 0.05% sodium azide
Concentration:	1.0 mg/ml
Purification:	Affinity chromatography (> 95% (by SDS-PAGE)
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	~ 130 kDa
Gene Name:	BUB1 mitotic checkpoint serine/threonine kinase B
Database Link:	Entrez Gene 12236 Mouse Entrez Gene 701 Human O60566



[View online »](#)

Background:

The N-end rule relates the in vivo half-life of a protein to the composition of its N-terminal residues. The N-end rule pathway is part of the ubiquitin system, which involves a three-step mechanism. Proteins targeted for degradation are bound on their N-terminal residue by Ubr1 (also designated E3 α and N-recognin), which catalyzes the covalent attachment of ubiquitin to the protein substrate. Two zinc finger domains and the RING-H2 finger domain of Ubr1 are essential for substrate recognition. Ubr1 is located on mouse chromosome 2 and on human chromosome 15 in the syntenic region. Ubr1 is a 220 kDa protein that is ubiquitously expressed in adult mouse, with the highest expression detected in skeletal muscle and heart. In mouse embryo, Ubr1 is primarily expressed in the branchial arches and in the tail and limb buds.

Synonyms:


Protein SSK1, BUBR1, MAD3L, BUB1 beta

Protein Families:

Druggable Genome, Protein Kinase, Stem cell - Pluripotency

Protein Pathways:

Cell cycle

Product images:


Western blot (WB) analysis of ubR1 antibody (Cat.-No.: [AP20400PU-N]) in extracts from HeLa cells treated with H₂O₂ 100uM 30.