

Product datasheet for **AP20409PU-N**

SnoN (SKIL) Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	Western blot: 1/500-1/1000. Immunohistochemistry on paraffin sections 1/500-1/1000.
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	This antibody detects endogenous levels of SnoN protein. (region surrounding Asp652)
Formulation:	Phosphate buffered saline (PBS), pH 7.2 State: Aff - Purified State: Liquid purified Ig fraction Preservative: 0.05% sodium azide
Concentration:	1.0 mg/ml
Purification:	Affinity-chromatography using epitope-specific immunogen; purity is > 95% (by SDS-PAGE)
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	~ 77 kDa
Gene Name:	SKI-like proto-oncogene
Database Link:	Entrez Gene 6498 Human P12757
Background:	SnoN may have a regulatory role in cell division or differentiation in response to extracellular signals. SnoN acts as a positive mediator of TGF beta induced transcription and cell cycle arrest in lung epithelial cells. The ability of SnoN to repress the growth inhibitory function of the Smad proteins is required for their transforming activity.



[View online »](#)

Synonyms: Ski-like protein, Ski-related protein, Ski-related novel protein N, Ski-related oncogene, SNO, SnoA, Snol

Product images:



Western blot analysis of SnoN antibody (Cat.-No.: AP20409PU-N) in extracts from HUVEC cells.