

Product datasheet for **AP20554PU-N**

Exonuclease 1 (EXO1) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IF, WB
Recommended Dilution:	Western blot: 1/500-1/1000. Immunofluorescence: 1/50-1/200.
Reactivity:	Human, Mouse
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	This antibody detects endogenous levels of Exo1 rotein. (region surrounding Lys86)
Formulation:	Phosphate buffered saline (PBS), pH~7.2 State: Aff - Purified State: Liquid purified Ig fraction (> 95% pure by SDS-PAGE). Preservative: 0.05% Sodium Azide
Concentration:	1.0 mg/ml
Purification:	Affinity Chromatography using epitope-specific immunogen.
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	~86,88 kDa
Gene Name:	exonuclease 1
Database Link:	Entrez Gene 26909 Mouse Entrez Gene 9156 Human Q9UQ84



[View online »](#)

Background:

Comparative evaluation of the expression patterns of the human and mouse genes, combined with previous biochemical and yeast genetic studies, indicate that the Exo1 (Exonuclease I) proteins are important contributors to chromosome processing during mammalian DNA repair and recombination. In mice, the mExo1 gene maps to distal chromosome 1, consistent with the recent mapping of the orthologous human HEX1/ hEXO1 gene to chromosome 1q42-q43. mExo1 is expressed prominently in testis, an area of active homologous recombination, and spleen, a prominent lymphoid tissue. In both mammalian and yeast systems, Exo1 is a 5'-3' double stranded DNA exonuclease that has previously been implicated in DNA mismatch repair (MMR). The mismatch repair (MMR) system ensures genome integrity by removing mispaired and unpaired bases that originate during replication. In humans, Exo1 interacts with MSH2 and MLH1 and has been proposed to be a redundant exonuclease in MM. In both mammalian and yeast systems, Exo1 plays a structural role in MMR and stabilizes multiprotein complexes containing a number of MMR proteins.

Synonyms:


Exonuclease I, EXOI, HEX1

Protein Families:

Druggable Genome, Stem cell - Pluripotency

Protein Pathways:

Mismatch repair

Product images:

Raw264.7 whole cell lysate
Exo1 (K86) pAb at 1:500 dilution

Western blot (WB) analysis of Exo1 antibody (Cat.-No.: AP20554PU-N) in extracts from raw264.7 cells.