


Product datasheet for **AP20644PU-N**

Human IgA Rabbit Polyclonal Antibody

Product data:


Product Type:	Secondary Antibodies
Product Name:	Human IgA Rabbit Polyclonal Antibody
Applications:	IHC, WB
Recommended Dilution:	Western blot: 1/500-1/1000. Immunohistochemistry on Paraffin Sections: 1/50-1/200.
Reactivity:	Human
Host:	Rabbit
Formulation:	PBS, pH~7.2 State: Aff - Purified State: Liquid purified Ig fraction (> 95% pure by SDS-PAGE). Preservative: 0.05% Sodium Azide
Concentration:	1.0 mg/ml
Purification:	Affinity Chromatography using epitope-specific immunogen.
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.

Product images:



Immunohistochemistry analyzes of IgA antibody (Cat.-No.: AP20644PU-N) in paraffin-embedded human lung carcinoma tissue.





Western blot analysis of IgA antibody (Cat.-No.: AP20644PU-N) in extracts from HeLa cells.