

Product datasheet for AP20652PU-N

ROCK1 Rabbit Polyclonal Antibody

Product data:

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	Western blot: 1/500-1/1000. Immunohistochemistry on paraffin sections: 1:50 - 1:200.
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	This antibody detects endogenous levels of Rock-1 protein. (region surrounding asn311)
Formulation:	Phosphate buffered saline (PBS), pH 7.2 State: Aff - Purified State: Liquid purified lg fraction Preservative: 0.05% sodium azide
Concentration:	1.0 mg/ml
Purification:	Affinity-chromatography using epitope-specific immunogen; purity is > 95% (by SDS-PAGE)
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	~ 158 kDa
Gene Name:	Rho associated coiled-coil containing protein kinase 1
Database Link:	<u>Entrez Gene 19877 MouseEntrez Gene 81762 RatEntrez Gene 6093 Human</u> <u>Q13464</u>




This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2024 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

GRIGENE ROCK1 Rabbit Polyclonal Antibody – AP20652PU-N


Background:	Rho, the Ras-related small GTPase, is responsible for the regulation of Actinbased cytoskeletal
	structures including stress fibers, focal adhesions and the contractile RING apparatus. Rho
	proteins function as molecular switches that are able to turn cytokinesis on and off. Although
	little is known about signaling downstream of Rho, a host of putative Rho effector proteins
	have been described, including rhophilin, Rhotekin, citron and the serine/threonine kinase, protein kinase N. Two additional Rho-activated serine/threonine kinases have been described, designated Rock-1 and Rock-2 (also referred to as Roka) for Rho-associated coil- containing protein kinase. Rock-1 and Rock-2 share a structural similarity with myotonic dystrophy kinase.
Synonyms:	ROCK-1, p160ROCK, NY-REN-35
Protein Families:	Druggable Genome, Protein Kinase
Protein Pathways:	Axon guidance, Chemokine signaling pathway, Focal adhesion, Leukocyte transendothelial migration, Pathogenic Escherichia coli infection, Regulation of actin cytoskeleton, TGF-beta

signaling pathway, Vascular smooth muscle contraction, Wnt signaling pathway

Product images:



Western blot (WB) analysis of Rock-1 antibody (Cat.-No.: AP20652PU-N) in extracts from COS7 cells.



Immunohistochemistry analyzes of Rock-1 antibody (Cat.-No.: AP20652PU-N) in paraffinembedded human brain tissue.

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2024 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US