

Product datasheet for **AP20937PU-N**

Estrogen Receptor 1 (ESR1) pSer104 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	Immunohistochemistry: 1/50 - 1/200.
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	p-ER-a (pSer104) pAb detects endogenous levels of p-ER-a protein
Formulation:	Phosphate buffered saline (PBS) with 15 mM sodium azide, approx. pH 7.2. State: Aff - Purified State: Liquid purified Ig fraction
Concentration:	1,0 mg/ml
Purification:	Affinity chromatography (> 95% (by SDS-PAGE)
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	estrogen receptor 1
Database Link:	Entrez Gene 13982 Mouse Entrez Gene 24890 Rat Entrez Gene 2099 Human P03372
Background:	The estrogen receptor-alpha (ER) is a broadly expressed transcription factor which controls a number of genes involved in cellular differentiation and proliferation. ER contains two transcriptionally active domains: the N-terminal domain (AF1) is constitutively active, while the C-terminal domain (AF-2) is hormone-dependent. ER is subject to posttranslational modification, including phosphorylation of serine residues 104, 106, 118, and 167. This antibody specifically recognizes human estrogen receptor-alpha when phosphorylated at serine 167, a site which influences AF-1 dependent transcriptional activity. This phosphorylation is catalyzed by the MAPK pathway, potentially through the 90 kDa ribosomal S6 kinase (pp90rsk1). Serine 167 is also a consensus Akt substrate.




[View online »](#)

Synonyms: ER alpha, Estradiol receptor, ESR1, ESR, NR3A1

Protein Families: Druggable Genome, Nuclear Hormone Receptor, Transcription Factors

Product images:



Immunohistochemistry (IHC) analyzes of p-ER-a (pSer104) pAb in paraffin-embedded human breast carcinoma tissue.