

## **Product datasheet for AP21252BT-N**

### **OriGene Technologies, Inc.** 9620 Medical Center Drive, Ste 200

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

# LOC100736962 Rabbit Polyclonal Antibody

#### **Product data:**

Product Type: Primary Antibodies

Applications: ELISA, ID, IF, IP, R, WB

**Recommended Dilution:** This product is intended for use in precipitating and non-precipitating antibody-binding

assays (such as e.g., ELISA and Western blotting and Immunofluorescence or Histochemical techniques), to prepare an insoluble immuno-affinity adsorbent, for labelling with a marker

of

the customer's own choice.

Working dilutions in non-precipitating antibody-binding techniques: 1/100-1/1,000.

Reactivity: Porcine
Host: Rabbit
Isotype: IgG

**Clonality:** Polyclonal

**Immunogen:** Esterase is isolated and purified from Porcine liver.

Freund's complete adjuvant is used in the first step of the immunization procedure.

**Specificity:** The antibody recognizes Esterase from Porcine liver.

The reagents were evaluated for potency, purity and specificity using most or all of the following techniques: Immunoelectrophoresis, Cross-Immunoelectrophoresis, single Radial

Immunodiffusion (Ouchterlony), block titration, ELISA, Immunoblotting and Enzyme

Inhibition.

Cross-reactivities against enzymes of other sources may occur but have not been

determined.

**Formulation:** PBS, pH 7.2 without preservatives and foreign proteins

Label: Biotin

State: Lyophilized hyperimmune IgG fraction

Molar radio: Biotin/IgG ~7.6

**Reconstitution Method:** Restore by adding 1.0 ml of sterile distilled water

**Concentration:** lot specific

**Purification:** Ammonium Sulphate Precipitation and Ion Exchange Chromatography

**Conjugation:** Biotin





#### LOC100736962 Rabbit Polyclonal Antibody - AP21252BT-N

Storage: Store the antibody lyophilized at 2-8°C and reconstituted at 2-8°C for one week or (in aliquots)

at -20°C for longer.

If a slight precipitation occurs upon storage, this should be removed by centrifugation.

**Stability:** Shelf life: one year from despatch.

Database Link: Entrez Gene 397478 Pig

Q29550

Synonyms: CES2, SES1, Liver carboxylesterase 1, HMSE, Egasyn, hBr1