

## Product datasheet for **AP21254BT-N**

### PFKM Sheep Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, ID, IF, IP, R, WB
Recommended Dilution:	This product is intended for use in precipitating and non-precipitating antibody-binding assays (such as e.g., ELISA and Western blotting and Immunofluorescence or Histochemical techniques). Working dilutions in non-precipitating antibody-binding techniques: 1/100-1/500.
Reactivity:	Rabbit
Host:	Sheep
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Fructose-6-Phosphate Kinase is isolated and purified from Rabbit muscle. Freund's complete adjuvant is used in the first step of the immunization procedure.
Specificity:	The antibody recognizes Fructose-6-Phosphate Kinase from Rabbit muscle. The reagents were evaluated for potency, purity and specificity using most or all of the following techniques: Immunoelectrophoresis, Cross-Immunoelectrophoresis, single Radial Immunodiffusion (Ouchterlony), block titration, ELISA, Immunoblotting and Enzyme Inhibition. Cross-reactivities against enzymes of other sources may occur but have not been determined.
Formulation:	PBS, pH 7.2 without preservatives and foreign proteins Label: Biotin State: Lyophilized hyperimmune IgG fraction Molar ratio: Biotin/IgG ~6.5
Reconstitution Method:	Restore by adding 1.0 ml of sterile distilled water
Concentration:	lot specific
Purification:	Ammonium Sulphate Precipitation and Ion Exchange Chromatography
Conjugation:	Biotin



[View online »](#)

<b>Storage:</b>	Store the antibody lyophilized at 2-8°C and reconstituted at 2-8°C for one week or (in aliquots) at -20°C for longer. If a slight precipitation occurs upon storage, this should be removed by centrifugation.
<b>Stability:</b>	Shelf life: one year from despatch.
<b>Database Link:</b>	<a href="#">P00511</a>
<b>Synonyms:</b>	Phosphofructokinase-M, Phosphohexokinase, Phosphofructokinase 1, PFK1, PFK-1, PFK-A, PFKX, GSD7