

## Product datasheet for **AP21368AF-N**

### PGM1 Sheep Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, ID, IF, IP, R, WB
Recommended Dilution:	Can be used in precipitating and non-precipitating antibody-binding assays (such as e.g., ELISA and Western blotting and immunofluorescence or histochemical techniques); to prepare an insoluble immuno-affinity adsorbent; for labelling with a marker of the customer's own choice. <u>Recommended working dilutions:</u> vary widely, but may be between 1/10 and 1/100.
Reactivity:	Rabbit
Host:	Sheep
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Phosphoglucomutase isolated and purified from rabbit muscle. Freund's complete adjuvant is used in the first step of the immunization procedure.
Specificity:	IgG fraction of polyclonal sheep antiserum to phosphoglucomutase from rabbit muscle. The reagents were evaluated for potency, purity and specificity using most or all of the following techniques: immunoelectrophoresis, cross-immunoelectrophoresis, single radial immunodiffusion (Ouchterlony), block titration, ELISA, immunoblotting and enzyme inhibition.
Formulation:	PBS, pH 7.2 without preservatives and foreign proteins State: Azide Free State: Lyophilized IgG fraction
Reconstitution Method:	Restore by adding 1.0 ml of sterile distilled water
Concentration:	lot specific
Purification:	Ammonium Sulphate Precipitation and Ion Exchange Chromatography
Conjugation:	Unconjugated
Storage:	Prior to and following reconstitution store the antibody at 2-8°C for one month or at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)

**Database Link:** [Entrez Gene 100009155 Rabbit P00949](#)

**Synonyms:** PGM 1, Phosphoglucomutase-1, Glucose phosphomutase 1