

## Product datasheet for **AP21372PU-N**

### PGAM2 Sheep Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, ID, IF, IP, R, WB
Recommended Dilution:	Can be used in precipitating and non-precipitating antibody-binding assays (such as e.g., ELISA and Western blotting and immunofluorescence or histochemical techniques); to prepare an insoluble immuno-affinity adsorbent; for labelling with a marker of the customer's own choice. <u>Recommended working dilutions:</u> vary widely, but may be between 1/100 and 1/2000.
Reactivity:	Rabbit
Host:	Sheep
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Phosphoglyceric mutase isolated and purified from rabbit muscle. Freund's complete adjuvant is used in the first step of the immunization procedure.
Specificity:	Affinity purified antibodies of polyclonal sheep antiserum to phosphoglyceric mutase from rabbit muscle. The reagents were evaluated for potency, purity and specificity using most or all of the following techniques: immunoelectrophoresis, cross-immunoelectrophoresis, single radial immunodiffusion (Ouchterlony), block titration, ELISA, immunoblotting and enzyme inhibition.
Formulation:	PBS, pH 7.2 stabilized with dextran without preservatives and foreign proteins State: Aff - Purified State: Lyophilized hyperimmune purified IgG fraction
Reconstitution Method:	Restore by adding 0.5 ml of sterile distilled water
Concentration:	lot specific
Purification:	Solid Phase Affinity Chromatography
Conjugation:	Unconjugated
Storage:	Prior to and following reconstitution store the antibody at 2-8°C for one month or at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)

**Gene Name:** phosphoglycerate mutase 2

**Database Link:** [Entrez Gene 5224 Human P15259](#)

**Synonyms:** PGAM 2, Phosphoglycerate mutase 2, PGAMM, PGAM-M

**Protein Families:** Druggable Genome

**Protein Pathways:** Glycolysis / Gluconeogenesis, Metabolic pathways