

Product datasheet for **AP21380BT-N**

PHKA1 Sheep Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, ID, IF, IP, R, WB
Recommended Dilution:	This product is intended for use in precipitating and non-precipitating antibody-binding assays (such as e.g., ELISA and Western blotting and immunofluorescence or histochemical techniques). <u>Working dilutions:</u> vary widely, but may be between 1/200 and 1/2000.
Reactivity:	Rabbit
Host:	Sheep
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Phosphorylase a isolated and purified from rabbit muscle. Freund's complete adjuvant is used in the first step of the immunization procedure.
Specificity:	Biotin-conjugated IgG fraction of polyclonal sheep antiserum to phosphorylase a from rabbit muscle. The reagents were evaluated for potency, purity and specificity using most or all of the following techniques: immunoelectrophoresis, cross-immunoelectrophoresis, single radial immunodiffusion (Ouchterlony), block titration, ELISA, immunoblotting and enzyme inhibition.
Formulation:	PBS, pH 7.2 without preservatives and foreign proteins Label: Biotin State: Lyophilized Hyperimmune IgG fraction Molar ratio: Biotin/ IgG ~6.0
Reconstitution Method:	Restore by adding 1.0 ml of sterile distilled water
Concentration:	lot specific
Purification:	Ammonium Sulphate Precipitation and Ion Exchange Chromat
Conjugation:	Biotin
Storage:	Prior to and following reconstitution store the antibody at 2-8°C for one month or at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)

Database Link: [Entrez Gene 100303771 Rabbit P18688](#)

Synonyms: MGC132604; PHKA