

## Product datasheet for **AP21401BT-N**

### Acetylcholinesterase Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, ID, IF, IP, R, WB
Recommended Dilution:	This product is intended for use in precipitating and non-precipitating antibody-binding assays (such as e.g., ELISA and Western blotting and immunofluorescence or histochemical techniques). <u>Recommended working dilutions:</u> vary widely, but may be between 1/200 and 1/2000.
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Acetylcholinesterase isolated and purified from electric eel. Freund's complete adjuvant is used in the first step of the immunization procedure.
Specificity:	The reagents were evaluated for potency, purity and specificity using most or all of the following techniques: immunoelectrophoresis, cross-immunoelectrophoresis, single radial immunodiffusion (Ouchterlony), block titration, ELISA, immunoblotting and enzyme inhibition.
Formulation:	PBS, pH 7.2 without preservatives and foreign proteins Label: Biotin State: Lyophilized IgG fraction Label: <u>Marker:</u> N-Hydroxysuccinimidobiotin. <u>Conjugation procedure:</u> A proprietary technique for the binding to biotin is used, followed by several purification steps. After each step activity and specificity are tested in a variety of techniques. The conjugate is lyophilised to assure stability and long shelf life Molar ratio: Biotin/ IgG ~6.1
Reconstitution Method:	Restore by adding 1.0 ml of sterile distilled water
Concentration:	lot specific
Purification:	Ammonium Sulphate Precipitation and Ion Exchange Chromatography
Conjugation:	Biotin
Storage:	Prior to and following reconstitution store the antibody at 2-8°C for one month or at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)

Stability:	Shelf life: one year from despatch.
Database Link:	<a href="#">O42275</a>
Synonyms:	AChE, Acetylcholine Esterase, Acetylcholine-Esterase