

Product datasheet for **AP21405BT-N**

Beta-galactose dehydrogenase Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, ID, IF, IP, R, WB
Recommended Dilution:	This antibody can be used in precipitating and non-precipitating antibody-binding assays (such as e.g., ELISA and Western blotting and Immunofluorescence or Histochemical techniques). <u>Recommended Dilutions:</u> Non-precipitating antibody-binding techniques: 1/100-1/2,000.
Reactivity:	<i>Pseudomonas fluorescens</i>
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Beta-Galactose Dehydrogenase is isolated and purified from <i>Pseudomonas fluorescens</i> . Freund's complete adjuvant is used in the first step of the immunization procedure.
Specificity:	This antibody is polyclonal immunologic reagents to beta-Galactose Dehydrogenase from <i>Pseudomonas fluorescens</i> . The reagents were evaluated for potency, purity and specificity using most or all of the following techniques: Immunoelectrophoresis, Cross-Immunoelectrophoresis, single Radial Immunodiffusion (Ouchterlony), block titration, ELISA, Immunoblotting and Enzyme Inhibition. Cross-reactivities against enzymes of other sources may occur but have not been determined.
Formulation:	PBS, pH 7.2 without preservatives and foreign proteins Label: Biotin State: Lyophilized hyperimmune IgG fraction Molar ratio: Biotin/IgG ~4.5
Reconstitution Method:	Restore by adding 1.0 ml of sterile distilled water
Concentration:	lot specific
Purification:	Ammonium Sulphate Precipitation and Ion Exchange Chromatography
Conjugation:	Biotin



[View online »](#)

Storage: Store the antibody lyophilized at 2-8°C and reconstituted at 2-8°C for one week or (in aliquots) at -20°C for longer.
If a slight precipitation occurs upon storage, this should be removed by centrifugation.

Stability: Shelf life: one year from despatch.