

## Product datasheet for **AP21469FC-N**

### Lactoferrin (LTF) Goat Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, ID, IF, IHC, IP
Recommended Dilution:	Can be used as reagent for the direct detection of lactoferrin in human cells, tissues and body fluids in immunofluorescence; as detection reagent in non-isotopic methodology and solid phase immunochemistry (e.g. ELISA). This immunoconjugate is not pre-diluted. The optimum working dilution of each conjugate should be established by titration before being used. Excess labelled antibody must be avoided because it may cause high unspecific background staining and interfere with the specific signal. <u>Working dilutions:</u> 1/20 - 1/80.
Reactivity:	Human
Host:	Goat
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen has been isolated from human milk. Freund's complete adjuvant is used in the first step of the immunization procedure.
Specificity:	In immunoelectrophoresis against human milk a single precipitin line is obtained. The antiserum does not react with any other protein component of human serum or plasma. <u>Cross-reactivity:</u> The antiserum does not cross-react with any other human plasma proteins as tested in gel-diffusion techniques. Inter-species cross-reactivity is a normal feature of antibodies to plasma or milk proteins, since homologous proteins of different species frequently share antigenic determinants. Cross-reactivity of this antiserum has not been tested in detail, however in double radial immunodiffusion a reaction with Rhesus monkey milk has been observed.
Formulation:	Label: FITC State: Lyophilized purified IgG fraction Label: Fluorescein Isothiocyanate Isomer 1 Absorption emission: 492 nm / 515 nm Molar ratio: Fluorochrome/IgG ~1.6
Reconstitution Method:	Restore by adding 1,0 ml of sterile distilled water
Concentration:	lot specific



[View online »](#)

<b>Purification:</b>	Hyperimmune antisera with strong precipitating activity are selected for Fractionation and purification of the IgG fraction containing the bulk of the defined antibody specificity.
<b>Conjugation:</b>	FITC
<b>Storage:</b>	Prior to and following reconstitution store the antibody at 2-8°C for one month or at -20°C for longer. Avoid repeated freezing and thawing.
<b>Stability:</b>	Shelf life: one year from despatch.
<b>Gene Name:</b>	lactotransferrin
<b>Database Link:</b>	<a href="#">Entrez Gene 4057 Human P02788</a>
<b>Background:</b>	Exocrine organs produce various secretions, each with its characteristic function. Proteins found in secretions may be divided into two groups: those specific for the particular secretion, and plasma proteins independent of the type of exocrine cells. Lactoferrin belongs to the first group. It is an iron containing protein with a molecular weight of 75,000 and it is antigenically different from transferrin. Lactoferrin has a slight antimicrobial action. Originally identified in milk, its presence has also been demonstrated in other secretions as saliva, semen and tears.
<b>Synonyms:</b>	LTF, LF, Lactoferrin, EC=3.4.21, Talalactoferrin
<b>Note:</b>	<u>Adsorption:</u> Immunoaffinity adsorbed using insolubilized fractions of human serum and lactoferrin-depleted human milk as required, to eliminate antibodies reacting with other human serum or milk proteins. The use of insolubilized adsorption antigens prevents the presence of excess adsorbent protein or immune complexes in the antiserum.
<b>Protein Families:</b>	Druggable Genome, Protease, Secreted Protein