

Product datasheet for AP21482AF-N

Mouse IgE (Fc specific) Goat Polyclonal Antibody

Product data:

Product Type: Secondary Antibodies Mouse IgE (Fc specific) Goat Polyclonal Antibody **Product Name: Applications:** ID, IF, IHC, IP, WB As unlabelled secondary antibody for indirect detection of IgE in mouse cell, tissue substrates Recommended Dilution: and body fluids in immunofluorescence and immunoenzyme assay methods; for the production of immunoconjugates with a selected marker; to prepare insoluble immunoaffinity adsorbents by coupling to an artificial carrier; as catching or detecting antibody reagent in non-isotopic assay methodology (e.g. ELISA) to identify and measure IgE in mouse serum or other body fluid. When applied in any immunocytochemical or histochemical staining procedure or solid phase coupling technique, the optimum concentration of the IgG preparation should be established. Recommended working dilutions: Histochemical and Cytochemical: 1/100 - 1/250. ELISA and comparable non-precipitating antibody-binding assays: 1/500 - 1/5000. Performance testing: Isotype specificity and quantitative specific recognition ability are further evaluated at the high level of sensitivity by direct singe and simultaneous double staining of different types of mouse cells, using counterstaining with reference conjugates with a different marker. The immunoconjugate is further tested in ELISA-type assays. **Reactivity:** Mouse Host: Goat Immunogen: Purified homogenous IgE isolated from mouse serum. Freund's complete adjuvant is used in the first step of the immunization procedure. Isotype: lgG Formulation: PBS, pH 7.2 without preservatives and foreign proteins State: Azide Free State: Lyophilized IgG fraction **Reconstitution Method:** Restore by adding 1.0 ml of sterile distilled water **Concentration:** 10,0 mg/ml **Purification:** Ammonium Sulphate Precipitation and Ion Exchange Chromatography **Conjugation:** Unconjugated



/iew online »

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

	Mouse IgE (Fc specific) Goat Polyclonal Antibody – AP21482AF-N
Storage:	Prior to and following reconstitution store the antibody at 2-8°C for one month or at -20°C for longer. Avoid repeated freezing and thawing.
Note:	Adsorption: Immunoaffinity adsorbed using insolubilized antigens as required to eliminate antibody activity reacting with other components of the immunoglobulin system and to any other serum protein. Special attention is given to the elimination of antibodies to the common Fab portion of immunoglobulins. The use of insolubilized adsorption antigens prevents the presence of excess adsorbent protein or immune complexes in the antiserum.

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US