

Product datasheet for **AP21482SU-N**

Mouse IgE (Fc specific) Goat Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgE (Fc specific) Goat Polyclonal Antibody
Applications:	ELISA, ID, IF, IHC, IP, WB
Recommended Dilution:	Can be used in precipitating techniques as immunoelectrophoresis and radial immunodiffusion to identify the presence of IgE in mouse serum and other body fluids or to determine its concentration. To prepare an immunoabsorbent for the purification of mouse IgE from serum or plasma. This antiserum is not intended for use in non-precipitating antibody-binding or other highly sensitive assays. For such use labelled and unlabelled cytochemical and immunoassay grade reagents meeting the necessary additional specificity and performance requirements are available. <u>Antibody titre:</u> Precipitin titre 1/32 when tested against monoclonal IgE containing mouse serum in agarblock immunodiffusion titration.
Reactivity:	Mouse
Host:	Goat
Immunogen:	Highly purified monoclonal IgE isolated from mouse serum. Freund's complete adjuvant is used in the first step of the immunization procedure.
Formulation:	State: Serum State: Lyophilized (Delipidated and heat inactivated) stable whole antiserum - No preservative added
Reconstitution Method:	Restore by adding 1,0 ml sterile distilled water. Dilutions may be prepared by adding PBS, pH 7.2
Concentration:	Total protein and IgG concentration in the antiserum are comparable to those of pooled normal goat serum. No foreign proteins added.
Conjugation:	Unconjugated
Storage:	Prior to and following reconstitution store the antibody at 2-8°C for one month or at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)

Note: Adsorption: Immunoaffinity adsorbed using insolubilized antigens as required, to eliminate antibodies cross-reacting with other components of the immunoglobulin system or reacting with other serum proteins. Special attention is given to the elimination of antibodies to the common Ig/Fab portion. The use of insolubilized adsorption antigens prevents the presence of excess adsorbent protein or immune complexes in the antiserum.