

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for AP22653PU-N

Activin Receptor Type IA (ACVR1) (N-term) Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	Immunohistochemistry on Paraffin Sections: 5 μg/ml. Western Blot: 1 μg/ml.
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	ACVR1 antibody was raised against aCVR1C antibody was raised against a 15 amino acid peptide near the amino terminus of the human ACVR1C.
Specificity:	This antibody reacts to Activin A Receptor, Type IA (ACVR1).
Formulation:	PBS containing 0.02% sodium azide State: Purified State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Immunoaffinity Chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C.
Stability:	Shelf life: one year from despatch.
Gene Name:	activin A receptor type 1
Database Link:	<u>Entrez Gene 11477 MouseEntrez Gene 79558 RatEntrez Gene 90 Human</u> <u>Q04771</u>
Synonyms:	Activin receptor type IA, ACVRLK2, ALK2, SKR1, ACRV1A
Protein Families:	Druggable Genome, ES Cell Differentiation/IPS, Protein Kinase, Transmembrane
Protein Pathways:	Cytokine-cytokine receptor interaction, TGF-beta signaling pathway




This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2024 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

Product images:




Western blot analysis of ACVR1C in human placenta tissue lysate with ACVR1C antibody at (A) 1 and (B) 2 ug/ml.



Human Adrenal (formalin-fixed, paraffinembedded) stained with ACVR1 antibody at 5 ug/ml followed by biotinylated goat anti-rabbit IgG secondary antibody, alkaline phosphatasestreptavidin and chromogen.

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2024 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US



Human Brain, Cortex (formalin-fixed, paraffinembedded) stained with ACVR1 antibody at 5 ug/ml followed by biotinylated goat anti-rabbit IgG secondary antibody, alkaline phosphatasestreptavidin and chromogen.

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2024 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US