

Product datasheet for **AP22936PU-N**

TRIM5 alpha (TRIM5) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	Immunohistochemistry on Paraffin Sections: 5 µg/ml. Western Blot: 1 - 2 µg/ml.
Reactivity:	Human, Monkey
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	TRIM5 antibody was raised against synthetic peptide - KLH conjugated
Specificity:	This antibody reacts to synthetic peptide corresponding to amino acids near the mid-region of rhesus monkey Tripartite Motif-containing 5 (TRIM5).
Formulation:	PBS, 0.02% sodium azide State: Aff - Purified State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Immunoaffinity Chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C to -70°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	tripartite motif containing 5
Database Link:	Entrez Gene 85363 Human Q9C035

[View online »](#)

Background:


TRIM5 is a member of a broad family of otherwise unrelated proteins defined by the presence of a tripartite motif containing a RING domain, a B-box type 1, and a B-box type 2, followed by a coiled-coil region. TRIM5 has six alternately spliced isoforms, the longest of which is the a variant which also contains a carboxy-terminal B30.2 (SPRY) domain (1). Expression of TRIM5a variants from humans, rhesus monkeys, and African green monkeys enabled resistance to infection by various retroviruses including HIV-1 (2,3), albeit at differing efficiencies. All TRIM5a variants could inhibit at least two different retroviruses, but not from those viruses isolated from the same species, suggesting that TRIM5a acts as a natural barrier to cross-species retrovirus transmission (3). While a function has not yet been assigned to TRIM5g, it is known that expression of TRIM5a variants from humans, rhesus monkeys, and African green monkeys enable resistance to infection by various retroviruses including HIV-1, albeit at differing efficiencies. Furthermore, the TRIM5d isoform appears to serve as a scaffold for the assembly of endogenous BTBD1 and BTBD2 proteins and also exhibits autoubiquitination activity in a RING finger- and UbCH5B-dependent manner. The TRIM5d isoform has been shown to colocalize with the topoisomerase I-interacting proteins BTBD1 and BTBD2 in punctate or elongated cytoplasmic bodies in several mouse and human cells where it appears to serve as a scaffold for the assembly of endogenous BTBD proteins. TRIM5d also exhibits autoubiquitination activity in a RING finger- and UbCH5B-dependent manner.

Synonyms:


RING finger protein 88

Protein Families:


Druggable Genome

Product images:


Western blot analysis of TRIM5a expression in human stomach (A), thymus (B), and uterus (C) cell lysate with TRIM5a antibody (AP22936PU-N) at 2 ug /ml.



IHC Image for AP22936PU-N Human Skeletal Muscle: Formalin-Fixed, Paraffin-Embedded (FFPE)



IHC Image for AP22936PU-N Human Prostate, stroma: Formalin-Fixed, Paraffin-Embedded (FFPE)