

## Product datasheet for **AP23243PU-N**

### Cyclin D2 (CCND2) (C-term) Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	Western blot: 1 µg/ml.
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide corresponding to a sequence at the C-terminal of human Cyclin D2
Specificity:	This antibody detects Cyclin D2 at C-term.
Formulation:	5mg BSA, 0.9mg NaCl, 0.2mg Na <sub>2</sub> HPO <sub>4</sub> , 0.05mg Thimerosal, 0.05mg Na <sub>3</sub> State: Aff - Purified State: Lyophilized Ig fraction
Reconstitution Method:	0.2ml of distilled water will yield a concentration of 500µg/ml.
Purification:	Immunogen affinity purified
Conjugation:	Unconjugated
Storage:	Store at 2 - 8 °C for up to one month or (in aliquots) at -20 °C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	cyclin D2
Database Link:	<a href="#">Entrez Gene 12444 Mouse</a> <a href="#">Entrez Gene 64033 Rat</a> <a href="#">Entrez Gene 894 Human P30279</a>



[View online »](#)

<b>Background:</b>	<p>Cyclin D2, also known as CCND2, is a human gene. The protein encoded by this gene belongs to the highly conserved cyclin family, whose members are characterized by a dramatic periodicity in protein abundance through the cell cycle. Cyclins function as regulators of CDK kinases. Different cyclins exhibit distinct expression and degradation patterns which contribute to the temporal coordination of each mitotic event. Cyclin D1, Cyclin D2 and Cyclin D3 are the members of the cyclin family. Cyclin D2 mapped to 12p131, since the CCND1 gene is on 11q13, this may be another bit of evidence of the homology of chromosomes 11 and 12. Choi D et al proved the expression of pseudogene cyclin D2 mRNA in the human ovary increases with age, which may be a novel marker for decreased ovarian function associated with the aging process.<sup>2</sup> And knockout studies of the homologous gene in mouse suggest the essential roles of this gene in ovarian granulosa and germ cell proliferation. High level expression of this gene was observed in ovarian and testicular tumors. Inaba et al. (1992) used murine cDNA clones for 3 cyclin D genes that are normally expressed during the G1 phase of the cell cycle to clone the cognate human genes. D-type cyclins (cyclins D1, D2, and D3) are regarded as essential links between cell environment and the core cell cycle machinery. Cyclin D2(CCND2) gene is assigned to 12p13.</p>
<b>Synonyms:</b>	Cyclin-D2, CCND2
<b>Protein Families:</b>	Druggable Genome
<b>Protein Pathways:</b>	Cell cycle, Focal adhesion, Jak-STAT signaling pathway, p53 signaling pathway, Wnt signaling pathway