

Product datasheet for **AP23453SU-S**

Domoic Acid Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA: 1/2,000.
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Immunization with Domoic Acid coupled to a carrier peptide (FID33) and additional immunization with Domoic Acid coupled to KLH
Specificity:	Domoic Acid.
Formulation:	State: Serum State: Lyophilized Serum
Reconstitution Method:	Restore in aqua bidest to initial volume.
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2-8°C for 6 months or at -20°C long term. After reconstitution store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C long term. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)

Product images:



Domoic Acid