

Product datasheet for **AP23827PU-N**

KIAA1881 (PLIN4) Goat Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC
Recommended Dilution:	ELISA: 1/2000. Immunohistochemistry on Paraffin Sections: 3.75 µg/ml.
Reactivity:	Human
Host:	Goat
Clonality:	Polyclonal
Immunogen:	Synthetic peptide (C-RTTDNGGEQTALSPQE) from an internal region of human KIAA1881
Specificity:	This antibody reacts to the internal region of PLIN4.
Formulation:	Tris saline buffer, pH 7.3 containing 0.5% BSA as stabilizer and 0.02% sodium azide as preservative State: Purified State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Immunoaffinity Chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	perilipin 4
Database Link:	Entrez Gene 729359 Human Q96Q06
Synonyms:	KIAA1881, Adipocyte protein S3-12, PLIN4



[View online »](#)

Product images:



Human Adipocytes: Formalin-Fixed, Paraffin-Embedded (FFPE)