

Product datasheet for **AP23847PU-N**

G RTP1 (N-term) Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	Immunohistochemistry on Paraffin Sections: 5 µg/ml. Western Blot: 1 - 2 µg/ml.
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide, 18 amino acid peptide from near the amino terminus of human GRTP1
Specificity:	This antibody reacts to the N-term of growth hormone regulated TBC protein 1.
Formulation:	PBS containing 0.02% sodium azide as preservative State: Purified State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Immunoaffinity Chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C.
Stability:	Shelf life: one year from despatch.
Gene Name:	growth hormone regulated TBC protein 1
Database Link:	Entrez Gene 66790 Mouse Entrez Gene 361180 Rat Entrez Gene 79774 Human Q5TC63
Synonyms:	Growth hormone-regulated TBC protein 1




[View online »](#)

Product images:



Western blot analysis of GTRP1 in K562 cell lysate with GTRP1 antibody at (A) 1 and (B) 2 ug/ml.



Human Kidney: Formalin-Fixed, Paraffin-Embedded (FFPE)