

Product datasheet for **AP23946PU-N**

MYBPC2 (C-term) Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	ELISA. Immunohistochemistry on Paraffin Sections: 5 µg/ml. Western Blot: 1 - 2 µg/ml.
Reactivity:	Human, Mouse
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	MYBPC2 antibody was raised against synthetic peptide - KLH conjugated
Specificity:	This antibody reacts to the C-term of Fast MyBP-C (MYBPC2).
Formulation:	PBS containing 0.02% sodium azide as preservative State: Purified State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Immunoaffinity Chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C.
Stability:	Shelf life: one year from despatch.
Gene Name:	myosin binding protein C, fast type
Database Link:	Entrez Gene 233199 Mouse Entrez Gene 4606 Human Q14324
Synonyms:	MYBPCF, Myosin-binding protein C



[View online »](#)

Product images:



Human Skeletal Muscle: Formalin-Fixed, Paraffin-Embedded (FFPE)