

Product datasheet for **AP26058PU-S**

Fbxl4 (1-5) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	Western blot: 1:500 - 1:1000.
Reactivity:	Drosophila
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Peptide sequence around aa.1~5 derived from Fruit fly F box and leucine-rich-repeat gene 4.
Specificity:	This antibody detects endogenous level of total F box and leucine-rich-repeat gene 4 protein.
Formulation:	Phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol State: Aff - Purified State: Liquid Ig fraction
Concentration:	lot specific
Purification:	Affinity-chromatography using epitope-specific peptide
Conjugation:	Unconjugated
Storage:	Store (in aliquots) at -20 °C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Database Link:	Q9VY46
Background:	A member of the family of mammalian F-box proteins. F-box proteins are an expanding family of eukaryotic proteins characterized by an approximately 40 aminoacid motif, the F box (so named because cyclin F was one of the first proteins in which this motif was identified) . Some F-box proteins have been shown to be critical for the controlled degradation of cellular regulatory proteins.
Synonyms:	FBL5, FBLX4
Note:	Molecular weight: 76 kDa



[View online »](#)