

Product datasheet for **AP26427AF-L**

Kcnd3 (55-64) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, WB
Recommended Dilution:	ELISA: 1:500,000 for detection of the peptide coated on ELISA plate. Western blot: 1:5,000 to detect the 74kDa mouse Kv4.3 in mouse brain tissue extract.
Reactivity:	Mouse
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	A synthetic peptide corresponding to amino acid 55-64 of mouse voltage-gated potassium channel Kv4.3, conjugated with KLH
Specificity:	This antibody is reactive to the 74kDa Kv4.3 in mouse brain tissue extract and other species with consensus Kv4.3 sequence (RTTLERYPDT).
Formulation:	0.01M PBS, pH 7.4 State: Azide Free State: Lyophilized Ig fraction
Reconstitution Method:	Restore with double distilled water to adjust the final concentration to 1.00mg/ml.
Purification:	Protein G affinity
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2-8°C for 6 month or at -20°C long term. After reconstitution store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C long term. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	potassium voltage-gated channel, Shal-related family, member 3
Database Link:	Entrez Gene 56543 Mouse Q9Z0V1



[View online »](#)

Background:

Kv4.3 is a member of the voltage-gated potassium channel family. Kv4.3 forms complex with K channel interacting protein and regulates the transient potassium efflux and repolarization in cardiac myocytes. Kv4.3 is highly expressed in brain tissue. It is a member of the rapid inactivating A-type Kv channel subfamily and is encoded by KCND3 gene in human. There are two isoforms of this protein with different sizes resulted from alternative splicing.

Synonyms:

Potassium voltage-gated channel subfamily D member 3, Voltage-gated potassium channel subunit Kv4.3

Product images:


Western blot analysis of extracts from mouse brain using anti-KV4.3 rabbit whole antiserum.