

## Product datasheet for AP31079HR-N

## 9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com

CN: techsupport@origene.cn

OriGene Technologies, Inc.

## **Transferrin Rabbit Polyclonal Antibody**

**Product data:** 

**Product Type:** Primary Antibodies

**Applications:** ELISA, ID, IF, IHC, IP, WB

**Recommended Dilution:** Can be used in enzyme-immunocytochemical and immunohistochemical techniques for the

detection of transferrin at the cellular and subcellular level in appropriately treated cell and

tissue substrates, as detection reagent in non-isotopic methodology and solid phase

immunochemistry (e.g. ELISA, Western blotting).

This immunoconjugate is not pre-diluted. The optimum working dilution of each conjugate should be established by titration before being used. Excess labelled antibody must be avoided because it may cause high unspecific background staining and interfere with the

specific signal.

**Recommended Working Dilutions:** 

Histochemical and Cytochemical: 1/100-1/500.

ELISA and comparable non-precipitating antibody-binding assays: 1/1,000-1/5,000.

Reactivity: Porcine

Host: Rabbit

Clonality: Polyclonal

**Immunogen:** Purified Transferrin isolated from pooled Swine serum.

Feund's complete adjuvant is used in the first step of the immunization procedure.

**Specificity:** The reactivity of the antiserum is restricted to Transferrin as tested in

Immunoelectrophoresis and radial immunodiffusion. A characteristic single precipitin line is obtained with normal serum which shows a reaction of identity with the precipitin line with

purified Transferrin.

**Cross-reactivity:** Inter-species cross-reactivity is a normal feature of antibodies to serum

proteins, since homologous proteins of differnt secies frequently share antigenic determinants. Cross-reactivity of this antiserum has not been tested in detail.

**Formulation:** PBS, pH 7.2 without preservatives and foreign proteins

Label: HRP

State: Lyphilized purified hyperimmune IgG fraction

Label: Horseradish Peroxidase Molar radio: Peroxidase/IgG ~1.7

**Reconstitution Method:** Restore with 1.0 ml distilled water.





## Transferrin Rabbit Polyclonal Antibody - AP31079HR-N

**Concentration:** lot specific

**Purification:** Hyperimmune antisera with strong precipitating activity are selected for fractionation and

purification of the IgG (7S) fraction containing the bulk of the defined antibody specificity. It is

free of other serum proteins as tested by immunoelectrophoresis.

Conjugation: HRP

Storage: Store lyophilized at 2-8°C and reconstituted at 2-8°C for one week or (in aliquots) at -20°C for

longer.

Avoid Repeated thawing and freezing.

**Stability:** Shelf life: one year from despatch.

Database Link: P09571

**Background:** Transferrin is a single polypeptide chain glycoprotein belonging to iron binding family of

proteins. It has a molecular weight of 81,5 kDa (porcine). It is synthesised in the liver and consists of two domains each having a high affinity reversible binding site for Fe3+. The iron is transported in blood and interstitial fluids to sites of use and disposal. Iron/transferrin is essential in haemoglobin synthesis and for certain types of cell division. Serum concentration rises in iron deficiency and pregnancy and falls in iron overload, infection and inflammatory

conditions.

**Synonyms:** Serotransferrin, Siderophilin

**Note:** Adsorption: Immunoaffinity adsorbed using insolubilized antigens as required, to eliminate

antibodies reacting with other serum proteins. The use of insolubilized adsorption antigens prevents the presence of excess adsorbent protein or immune complexes in the antiserum. **Conjugation procedure:** Conjugation is carried out using a proprietary modification of the periodate technique for the binding to peroxidase, followed by several purification steps. After each step activity and specificity are tested in a variety of techniques. The conjugate is

lyophilized to assure stability and long shelf life.