

## Product datasheet for **AP31195PU-N**

### **FOXO4 (N-term) Rabbit Polyclonal Antibody**


#### **Product data:**

Product Type:	Primary Antibodies
Applications:	ELISA, IF, IHC, WB
Recommended Dilution:	<b>ELISA.</b> <b>Immunocytochemistry.</b> <b>Immunohistochemistry on Paraffin Sections:</b> 10 µg/ml. <b>Western Blot:</b> 0.5 - 1 µg/ml.
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	FOXO4 antibody was raised against synthetic peptide - KLH conjugated
Specificity:	This antibody reacts to the N-term of FOXO4.
Formulation:	PBS containing 0.02% sodium azide as preservative State: Purified State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Immunoaffinity chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C.
Stability:	Shelf life: one year from despatch.
Gene Name:	forkhead box O4
Database Link:	<a href="#">Entrez Gene 54601 Mouse</a> <a href="#">Entrez Gene 302415 Rat</a> <a href="#">Entrez Gene 4303 Human</a> <a href="#">P98177</a>
Synonyms:	Forkhead box protein O4
Protein Families:	Transcription Factors




[View online »](#)

**Product images:**



Human Adrenal: Formalin-Fixed, Paraffin-Embedded (FFPE)



Human Vessels: Formalin-Fixed, Paraffin-Embedded (FFPE)