

## Product datasheet for **AP31201PU-N**

### **SYK (313-339) Rabbit Polyclonal Antibody**


#### **Product data:**

|                       |   |
|-----------------------|---|
| Product Type:         | Primary Antibodies  |
| Applications:         | IHC, IP, WB   |
| Recommended Dilution: | <b>Immunohistochemistry on Paraffin Sections:</b> 1/50.<br><b>Immunoprecipitation.</b><br><b>Western Blot.</b>                                  |
| Reactivity:           | Human   |
| Host:                 | Rabbit  |
| Isotype:              | IgG   |
| Clonality:            | Polyclonal  |
| Immunogen:            | Synthetic peptide derived from aa 313-339 of human Syk protein  |
| Specificity:          | This antibody reacts to SYK.  |
| Formulation:          | 10 mM PBS, pH 7.4 containing BSA as stabilizer and sodium azide as preservative<br>State: Purified<br>State: Liquid purified Ig fraction        |
| Purification:         | Immunoaffinity chromatography   |
| Conjugation:          | Unconjugated  |
| Storage:              | Store the antibody undiluted at 2-8°C.<br><b>DO NOT FREEZE!</b>   |
| Stability:            | Shelf life: one year from despatch.   |
| Gene Name:            | spleen tyrosine kinase  |
| Database Link:        | <a href="#">Entrez Gene 6850 Human P43405</a>   |
| Synonyms:             | Tyrosine-protein kinase SYK, Spleen tyrosine kinase   |
| Protein Families:     | Druggable Genome, Protein Kinase  |
| Protein Pathways:     | B cell receptor signaling pathway, Fc epsilon RI signaling pathway, Fc gamma R-mediated phagocytosis, Natural killer cell mediated cytotoxicity |



[View online »](#)

**Product images:**



Human Tonsil: Formalin-Fixed, Paraffin-Embedded (FFPE)