

Product datasheet for **AP31239PU-N**


BANP Rabbit Polyclonal Antibody

Product data:


| | |
|-----------------------|--|
| Product Type: | Primary Antibodies |
| Applications: | ELISA, IHC, WB |
| Recommended Dilution: | ELISA. Immunohistochemistry on Paraffin Sections: 5 µg/ml. Western Blot: 1 - 2 µg/ml. |
| Reactivity: | Human, Mouse, Rat |
| Host: | Rabbit |
| Clonality: | Polyclonal |
| Immunogen: | Synthetic peptide - KLH conjugated |
| Specificity: | This antibody reacts to the internal region of BANP. |
| Formulation: | PBS containing 0.02% sodium azide as preservative State: Purified State: Liquid purified Ig fraction |
| Concentration: | lot specific |
| Purification: | Immunoaffinity chromatography |
| Conjugation: | Unconjugated |
| Storage: | Store the antibody undiluted at 2-8°C. |
| Stability: | Shelf life: one year from despatch. |
| Gene Name: | BTG3 associated nuclear protein |
| Database Link: | Entrez Gene 53325 Mouse Entrez Gene 292064 Rat Entrez Gene 54971 Human Q8N9N5 |
| Synonyms: | BEND1 |




[View online »](#)

Product images:


Human Brain, Cerebellum: Formalin-Fixed, Paraffin-Embedded (FFPE)



Human Placenta: Formalin-Fixed, Paraffin-Embedded (FFPE)



Human Spleen: Formalin-Fixed, Paraffin-Embedded (FFPE)



Human Thymus: Formalin-Fixed, Paraffin-Embedded (FFPE)