

Product datasheet for AP31279PU-N

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

Apolipoprotein H (APOH) Goat Polyclonal Antibody

Product data:

Product Type: Primary Antibodies

Applications: ELISA, ID, IHC

Recommended Dilution: ELISA.

Immunodiffusion.

Immunoelectrophoresis.

Immunohistochemistry on Paraffin Sections: 5 µg/ml.

Reactivity: Canine, Human, Porcine, Rat

Host: Goat Isotype: IgG

Clonality: Polyclonal

Immunogen: Purified beta2-Glycoprotein-I (beta2GP-I) from human plasma. This protein is also known as

apolipoprotein-H.

Specificity: This antibody recognizes Human beta2GP-I as demonstrated by immunoelectrophoresis. A

single precipitin arc is observed with normal plasma. No arc was seen with beta2GP-I

deficient plasma.

Formulation: 10mM HEPES, pH 7.4, 0.15M sodium chloride, 50% glycerol.

State: Aff - Purified

State: Liquid purified Ig fraction

Concentration: lot specific

Purification: Ammonium Sulfate Precipitation

Conjugation: Unconjugated

Storage: Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer.

Avoid repeated freezing and thawing.

Stability: Shelf life: one year from despatch.

Gene Name: apolipoprotein H

Database Link: Entrez Gene 287774 RatEntrez Gene 350 Human

P02749






Apolipoprotein H (APOH) Goat Polyclonal Antibody - AP31279PU-N

Synonyms: Beta-2-glycoprotein I, Beta2GPI, B2GPI, ApoH, Apo-H, APC Inhibitor, Anticardiolipin Cofactor,

B2G1

Protein Families: Druggable Genome, Secreted Protein

Product images:



Human Liver: Formalin-Fixed, Paraffin-Embedded (FFPE)