

## Product datasheet for **AP31390AF-N**

### Vancomycin Sheep Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA.
Host:	Sheep
Clonality:	Polyclonal
Immunogen:	Vancomycin - KLH Conjugate
Specificity:	This antibody reacts to Vancomycin. Minimal cross-reactivity with related compounds.
Formulation:	State: Azide Free State: Liquid purified Ig fraction
Conjugation:	Unconjugated
Storage:	Store the antibody at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Note:	Titer: 1/7,000



[View online »](#)