

Product datasheet for **AP31470SU-N**

Equine Serum Proteins Rabbit Polyclonal Antibody

Product data:

Product Type: Primary Antibodies

Applications: ID, IP, R

Recommended Dilution: In precipitating techniques as immunoelectrophoresis and radial immunodiffusion (Ouchterlony) to identify the serum protein pattern, or the presence or absence of an individual component.

To evaluate the purity of an isolated serum protein including immunoglobulins. Since immunoprecipitation depends on a correct antigen/antibody concentration ratio (zone of equivalence) in the gel medium, the protein analysis by immunoelectrophoresis of serum or any other biological fluid or protein fraction should include different proportions of the reactants. It is not possible to obtain an optimal protein pattern in a single analysis. The electroendosmosis effect of different types of agar on proteins with a different net charge can be used to optimize the resolution power of the test system. Agar Nordic Nr. 2 contains sufficient positively charged ions to optimize the resolution of the proteins in the beta-gamma regions, while the alpha regions will become more dense. Highly purified Agar with low electroendosmosis favours the resolution of the proteins in the alpha regions, while the major components in the beta-gamma region can still be identified.

Directions for use In immunoelectrophoresis use 2 ml serum or equivalent against 120 ml antiserum. In double radial immunodiffusion (Ouchterlony) use a rosette arrangement with 10 ml antiserum in 3 mm diameter center well and 2 ml serum samples (neat and serially diluted in 2 mm diameter peripheral wells).

Reactivity: Equine

Host: Rabbit

Clonality: Polyclonal

Immunogen: Pooled whole horse serum and partly purified serum fractions.
Freund's complete adjuvant is used in the first step of the immunization procedure.



[View online »](#)

Specificity:	<p>In immunoelectrophoresis against pooled serum and enriched serum proteins fractions precipitation can be observed of not less than 12 individual proteins components.</p> <p>Cross-reactivity: Inter-species cross-reactivity is a normal feature of antibodies to animal proteins since homologous proteins of different species frequently share antigenic determinants.</p> <p>This antiserum has not been adsorbed for such cross-reactivity. Consequently it is not species-specific.</p>
Formulation:	<p>State: Serum</p> <p>State: Delipidated, heat inactivated, lyophilized, stable whole antiserum.</p> <p>Preservative: None</p>
Reconstitution Method:	<p>Restore the lyophilized antiserum by adding 1 ml sterile distilled water. Dilutions may be prepared by adding phosphate buffered saline (PBS, pH 7.2).</p>
Concentration:	<p>Total protein and IgG concentrations in the antiserum are comparable to those of pooled normal rabbit serum. No foreign proteins added.</p>
Conjugation:	<p>Unconjugated</p>
Storage:	<p>Store lyophilized at 2-8°C for 6 months or at -20°C long term.</p> <p>After reconstitution store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C long term.</p> <p>If a slight precipitation occurs upon storage, this should be removed by centrifugation. It will not affect the performance of the antiserum</p> <p>Avoid repeated freezing and thawing.</p>
Stability:	<p>Shelf life: one year from despatch.</p>
Synonyms:	<p>Horse Serum Proteins</p>