

## Product datasheet for AP31521SU-N

## **C3** Rabbit Polyclonal Antibody

**Product data:** 

**Product Type: Primary Antibodies** 

**Applications:** ID, IHC, IP

Recommended Dilution: Can be used in precipitating techniques as immunoelectrophoresis and single and double

> radial immunodiffusion (Mancini, Ouchterlony) to identify the presence of complement C3c or to determine its concentration. The presence of non-precipitating antibodies has not been assayed. This does not exclude the use of the antiserum in non-precipitating antibodybinding techniques if proper controls are included. Determinations of individual complement

OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com

components can be very useful in defining the exact location of a defect.

Also suitable for Immunohistochemistry on Paraffin Sections.

Directions for use.

In immunoelectrophoresis use 2 µl human plasma or equivalent against 120 µl antiserum. In double radial immunodiffusion use a rosette arrangement with 10 µl antiserum in 3 mm diameter center well and 2 µl plasma samples (neat and serially diluted) in 2 mm diameter

peripheral wells.

In single radial immunodiffusion use 1 percent antiserum in the gel.

Immunohistochemistry on Paraffin Sections: 1/200.

Reactivity: Human Host: Rabbit Clonality: Polyclonal

Immunogen: C3c is the major fragment resulting from the C3 cleavage by C3 convertase and factor i. It is

composed of an intact beta chain bound to two fragments of the alpha chain. The protein is

isolated and purified from pooled normal human serum by precipitation techniques,

followed by chromatographical methods.

Freund's complete adjuvant is used in the first step of the immunization procedure.





Formulation:

## C3 Rabbit Polyclonal Antibody - AP31521SU-N

**Specificity:** In immunoelectrophoresis against fresh human serum, a single precipitin line is obtained in

the beta-1 region representing native C3. Against serum containing partly activated C3, a

precipitin line is obtained which extends from the beta-1 into the alpha-2 region,

determinants. Crossreactivity of this antiserum has not been tested in detail.

demonstrating a gradient. In old serum containing totally activated C3 a single precipitin line in the alpha-2 region is obtained. Antisera to C3c can also react with the fragments C3b, C3bi and smaller fragments, since they all carry antigenic determinants of the C3c domain. The product does not react with any other proteins component of human serum or plasma. <a href="Cross-reactivity">Cross-reactivity:</a> Inter-species cross-reactivity is a normal feature of antibodies to plasma proteins, since homologous proteins of different species frequently share antigenic

No preservative added. - No foreign proteins added.

State: Serum

State: Delipidated, heat inactivated, lyophilized, stable whole serum.

**Reconstitution Method:** Restore with 1 ml sterile distilled water.

Conjugation: Unconjugated

**Storage:** Prior to reconstitution store at 2-8°C.

Following reconstitution store undiluted at 2-8°C for one week

or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.

**Stability:** Shelf life: one year from despatch.

Gene Name: complement component 3

Database Link: Entrez Gene 717 Human

P01026

**Synonyms:** CPAMD1, Complement component 3

**Note:** Adsorption: Immunoaffinity adsorbed using insolubilized antigens as required, to eliminate

antibodies reacting with other human serum proteins. The use of insolubilized adsorption antigens prevents the presence of excess adsorbent protein or immune complexes in the

antiserum.

**Protein Families:** Druggable Genome, Protease, Secreted Protein

**Protein Pathways:** Complement and coagulation cascades, Systemic lupus erythematosus