

## Product datasheet for **AP31649SU-N**

### Human IgG Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgG Rabbit Polyclonal Antibody
Applications:	ID, IP
Recommended Dilution:	<b>Immunoprecipitation.</b> <b>Immunodiffusion.</b> <u>For Ouchterlony:</u> Sample Size: 15 µl Distance between wells: 5 mm Test time: 4 days at 37°C Antigens tested: Normal human serum, purified human IgG, normal mouse serum. Results: Titre with NHS > 1:1000.
Reactivity:	Human
Host:	Rabbit
Formulation:	State: Serum State: Antiserum containing 0.05% Sodium Azide as preservative
Conjugation:	Unconjugated
Storage:	Store the antibody (in aliquots) at -20°C.



[View online »](#)