

Product datasheet for **AP31651SU-N**

PMN Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	CT
Recommended Dilution:	Cytotoxic Assays. In-vivo Depletion. Test Type: For Cytotoxic Antibodies: Modified Colormetric Microtiter Assay. (1) Results: Antisera of this lot (6326) diluted 1:100 exhibits >85% cytotoxicity on Mouse PMNs. For Agglutinating Antibodies: Antisera dilutions in RPMI-1640 incubated with target cells at 4-8°C for 1hr. Agglutination determined by microscopic observations. Results: Antisera of this lot (6326) strongly agglutinates Mouse PMNs at dilutions 1:200.
Reactivity:	Mouse
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	This antibody is an antiserum directed against Mouse PMNs (Polymorphonuclear Leukocytes).
Formulation:	State: Serum State: Liquid Antiserum. Preservative: None
Conjugation:	Unconjugated
Storage:	Upon receipt, store undiluted (in aliquots) at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Background:	Granulocytes are a category of white blood cells characterized by the presence of granules in their cytoplasm. They are also called polymorphonuclear leukocytes (PMN or PML) because of the varying shapes of the nucleus, which is usually lobed into three segments. Granulocytes or PMN are released from the bone marrow by the regulatory complement proteins.
Synonyms:	Polymorphonuclear Leukocytes



[View online »](#)