

Product datasheet for **AP31676PU-N**

C1QA Goat Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA.
Reactivity:	Human
Host:	Goat
Isotype:	IgG
Clonality:	Polyclonal
Specificity:	This antibody reacts to Complement C1q.
Formulation:	10 mM Sodium Phosphate, 0.15 M Sodium Chloride, pH 7.2 State: Aff - Purified State: Liquid 0.2µm filtered purified Ig fraction Preservative: 0.05% (w/v) Sodium Azide
Concentration:	lot specific
Purification:	Affinity purified using solid-phase human Complement C1q
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C.
Stability:	Shelf life: one year from despatch.
Gene Name:	complement component 1, q subcomponent, A chain
Database Link:	Entrez Gene 727 Human Entrez Gene 712 Human P02745
Synonyms:	Complement C1q, Complement 1q
Protein Families:	Druggable Genome
Protein Pathways:	Complement and coagulation cascades, Prion diseases, Systemic lupus erythematosus



[View online »](#)